



BIOTEQ®
Improving lives



**PTA BALLOON
DILATATION CATHETER**
(Percutaneous Transluminal Angioplasty)



BIOTEQUE CORPORATION




The BIOTEQ[®] Over the wire balloon dilatation catheter offers high performance femoral artery dilatation. For the treatment of obstructive lesions on native or synthetic arteriovenous dialysis fistula, this product delivers cost-effective solution for all patients.



Low compliant balloon with material reinforced Nylon provides powerful penetration to obstructive lesions.




Improved retrieving profile expedites easier crossing through stenosis and occlusions.






The BIOTEQ® Rapid Exchange balloon dilatation catheter is designed for applications in the peripheral vasculature, including iliac, femoral, ilio-femoral, popliteal, and renal arteries.



Semi-compliance nylon material balloon exhibits relatively high tensile strength, controlled compliance and with lower resistance.



With a side guidewire aperture design, balloon catheter can be operated by a single user when exchanging guidewire.



Superior triple balloon refolding characteristics.




Low profile offers excellent lesion crossibility.



Technical information – specifications

□ Rapid Exchange




Product Code	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Length (cm)	Recommended Guidewire Size (inch)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Recommended Introducer Sheath Size (Fr)
BT-BC-R4-3-15-135	3.0	15	135	0.014	6	16	4
BT-BC-R4-3-20-135	3.0	20	135	0.014	6	16	4
BT-BC-R8-4-40-40	4.0	40	40	0.018	6	18	4
BT-BC-R8-5-40-40	5.0	40	40	0.018	6	18	4
BT-BC-R8-6-40-40	6.0	40	40	0.018	6	18	4

□ Over the Wire




Dual luman type



Product Code	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Length (cm)	Recommended Guidewire Size (inch)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Recommended Introducer Sheath Size (Fr)
BT-BC-D5-5-40-75	5.0	40	75	0.035	6	16	6
BT-BC-D5-6-40-75	6.0	40	75	0.035	6	16	6
BT-BC-D5-7-40-75	7.0	40	75	0.035	6	16	7
BT-BC-D5-8-40-75	8.0	40	75	0.035	6	16	7
BT-BC-D5-5-40-120	5.0	40	120	0.035	6	16	6
BT-BC-D5-6-40-120	6.0	40	120	0.035	6	16	6
BT-BC-D5-7-40-120	7.0	40	120	0.035	6	16	7
BT-BC-D5-8-40-120	8.0	40	120	0.035	6	16	7

Coaxial Type



Product Code	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Length (cm)	Recommended Guidewire Size (inch)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Recommended Introducer Sheath Size (Fr)
BT-BC-C8-4-20-40	4.0	20	40	0.018	6	20	4
BT-BC-C8-4-40-40	4.0	40	40	0.018	6	20	4
BT-BC-C8-5-20-40	5.0	20	40	0.018	6	20	4
BT-BC-C8-5-40-40	5.0	40	40	0.018	6	20	4
BT-BC-C8-6-20-40	6.0	20	40	0.018	6	18	4
BT-BC-C8-6-40-40	6.0	40	40	0.018	6	18	4