



BIOTEQ®

Improving lives

LONG TERM PIGTAIL DRAINAGE CATHETER




- Designed for percutaneous drainage of abscesses applied in gastroenterological system, urological system and other applicable drainage abscesses.
- Do not use continuously for more than 90 days. The catheter should be evaluated by the physician on or before 90 days post-placement periodically. Made from a soft, biocompatible radiopaque polyurethane (TPU).
- If required the operator can choose different drainage sets according to the type of accumulated fluid and place of accumulation.


CATHETER TIP DESIGNATION




Straight type




Pigtail type




(V)=(J)
Curved type




(V)=(MN)
Mini-pigtail type



(V)=(M)
Malecot type



(V)=(CL)
Closed loop
pigtail type





(V)=(MCL)
Mini Closed loop
pigtail type

DESCRIPTION/ FEATURES

- Used for drainage of nephrostomy, abscess, and other multipurpose drainage.
- Without or without safety string lock.
- Seven different tip design to meet your procedural needs.
- Tapered tip for smooth introduction.
- Echo enhanced needle for increased opacity.

LONG TERM PIGTAIL CATHETER SET-DIRECT ACCESS



Product Code

Straight tip: BT-PDL0-XXYY

Pigtail tip: BT-PDL1-XXYY-W

Other tips: BT-PDL1-XXYY-W(V)

XX = 05*, 06, 07, 08, 09, 10, 12, 14, 16 FR.

W = With string lock

*Without string lock

YY = 20, 25, 30, 35, 40cm



(V) = See catheter tip designation

Contents:

Drainage catheter x 1 pc,

Puncture Needle x 1pc

LONG TERM PIGTAIL CATHETER SET-SEDLINGER



Product Code

Straight tip: BT-PDL2-XXYY

Pigtail tip: BT-PDLS-XXYY-W

Other tips: BT-PDLS-XXYY-W(V)

XX = 05*, 06, 07, 08, 09, 10, 12, 14, 16 FR.

YY = 20, 25, 30, 35, 40, 45, 50cm

W = With string lock

(V) = See catheter tip designation

*Without string lock



Contents:

Drainage catheter x 1 pc,

Flexible stiffening cannula x 1pc,

Metal stiffening cannula x 1pc

LONG TERM PIGTAIL CATHETER



Product Code

Straight tip: BT-PDL3-XXYY

Pigtail tip: BT-PDL4-XXYY-W

XX = 05*, 06, 07, 08, 09, 10, 12, 14, 16 FR.

YY = 20, 25, 30, 35, 40, 45, 50cm

W = With string lock

*Without string lock

Contents:

Drainage catheter x 1 pc

For more information, please contact:
sales@bioteq.com.tw

BIOTEQUE CORPORATION
5F-6, No. 23, Sec. 1, Chang-An E. Rd.,
Taipei 104, Taiwan.
Phone: (886)-2-2571-0269, Fax: (886)-2-2536-1967
<http://www.bioteq.com.tw>

ISO 13485, GMP, CE Certified
© Bioteque Corp. Mar. 2016
Brochure number: BT-PDL-01